


Compatibility mode for Windows 7

Follow the directions and pictures to run BioExplorer in Windows 7.

Right-click on the BioExplorer application icon and select Properties

Computer\Local disk (C:)\Program Files (x86)\BioExplorer > right-click BioExplorer > Properties

First go to Computer\Local disk (C:)\Program Files (x86)\BioExplorer and then in BioExplorer (see 1st image) right click on BioExplorer and scroll down to Properties at the bottom (see Properties highlighted in blue).



Tick "Run the program in compatibility mode for" and select Windows 7 from the dropdown menu.

Properties > Compatibility > Compatibility Mode > Run this Program in Compatibility mode for > Windows 7

Go to the Compatibility tab at the top and check the box "Run this Program in Compatibility mode for:" and choose from the drop down menu – Windows 7 See the image in the BioExplorer Properties box. Then check OK.

